

# Glycerol Kinase“Amano”2

(ATP : Glycerol 3-phosphotransferase, EC 2.7.1.30)

Glycerol Kinase“Amano”2 is a glycerol kinase preparation manufactured by a unique fermentation process using a selected strain of bacteria.

## Catalysis



## Specification and Preparation

Activity:	Glycerol kinase activity	$\geq 10$ u/mg (Amano method)
Contaminants:	ATPase activity	$\leq 1.0 \times 10^{-2}\%$
	Catalase activity	$\leq 1.0 \times 10^{-1}\%$
	Hexokinase activity	$\leq 5.0 \times 10^{-2}\%$
	NADH oxidase activity	$\leq 1.0 \times 10^{-2}\%$
Appearance:	White to light yellowish white powder, lyophilized	
Additive:	Sodium glutamate	

## Characteristics

1. Molecular weight: 82,000 (Gel filtration)
2. Isoelectric point: 4.3
3. Km:  $3.5 \times 10^{-5}$  M (Glycerol)  
 $3.7 \times 10^{-5}$  M (ATP)
4. Optimum pH: 10.0
5. pH stability: 5.0-9.5 (37°C, 1hr)
6. Optimum temperature: 50°C
7. Thermal stability: up to 50°C (pH 7.0, 15 min)
8. Inhibitors: SH reagents
9. Stabilizer: ATP
10. Application: Used for the enzymatic determination of triglyceride in serum by coupling with glycerophosphate dehydrogenase (G-3-PDH) and L- $\alpha$ -glycerophosphate oxidase (GPO) or pyruvate kinase and lactate dehydrogenase in clinical diagnosis.

## Expiration (Storage)

24 months from the date of analysis when stored at -20°C or below in a dry place under sealed conditions.

The information and recommendations contained herein are to the best of our knowledge reliable according to the current scientific and technical level. However, depending upon use method and/or condition, nothing herein is to be construed as a warranty or representation in respect otherwise, including freedom from patent infringement. Users shall make their own test and investigation for their particular purpose. We do not accept any liability for any loss, damage or infringement arising from the use of information and recommendations contained herein.


## Safe Handling

1. Do not inhale.
2. In case of direct contact with skin or eyes, immediately wash or rinse with plenty of water.
3. Please refer to SDS for more details.


## General properties

The following results demonstrate the activities of enzyme solution prepared in various buffers. Enzyme activity may vary under different experimental conditions.


### pH and Activity




### pH stability



### Temperature and Activity




### Thermostability



### Substrate Specificity


Substrate	Relative Activity (%)
Glycerol	100
Dihydroxyacetone	97
DL-Glyceraldehyde	55

### Stability (powder form)




# Assay method of Glycerol kinase activity

## Principle



\*G3PDH: Glycerol-3-phosphate dehydrogenase



The appearance of NADH is measured at 340 nm by spectrophotometry.

## Unit Definition

One unit is defined as the enzyme quantity which produces one  $\mu$  mole of NADH per minute under the conditions described below.

## Reagents

- A. 0.2M Glycine-hydrazine buffer (containing 2 mM  $\text{MgCl}_2$ , pH 9.8)  
Weigh 15 g of glycine, 51.1 g of hydrazine hydrate (98%, FUJIFILM Wako Pure Chemical Corporation) and 406 mg of  $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$  then dissolve in 800 ml of deionized water. Adjust the pH to 9.8 at 25°C with 1N KOH and fill up to 1000 ml with deionized water. (Can be used for 3 months if kept refrigerated)
- B. 0.02M Phosphate buffer ( $\text{KH}_2\text{PO}_4$ - $\text{Na}_2\text{HPO}_4$ , pH 7.0)
- C. NAD/ATP solution  
Weigh 10 mg of  $\beta$ -NAD (Oriental Yeast Co., Ltd.) and 30 mg of  $\text{ATP} \cdot 2\text{Na}$  (Oriental Yeast Co., Ltd.), then dissolve in 20 ml of 0.2M Glycine-hydrazine buffer (A). (Can be used for 5 days if kept refrigerated)
- D. Substrate solution  
Weigh 92 mg of glycerol (FUJIFILM Wako Pure Chemical Corporation) and dissolve in 0.2 M Glycine-hydrazine buffer (A). Fill up to 100 ml with 0.2 M Glycine-hydrazine buffer (A). (Can be used for 1 month if kept refrigerated)
- E. G3PDH solution  
(Produced by Roche Diagnostics GmbH)
- F. Diluent  
Weigh 303 mg of  $\text{ATP} \cdot 2\text{Na}$  (Oriental Yeast Co., Ltd.) and dissolve in 0.02M Phosphate buffer (B). Fill up to 500 ml with 0.02M Phosphate buffer (B). (Can be used for 10 days if kept refrigerated)
- G. Enzyme solution  
Weigh out Glycerol Kinase "Amano"2 and dissolve in chilled Diluent (F). Enzyme solution should be prepared so that the value of  $\Delta\text{OD}/\text{minute}$  becomes in the range of  $0.030 \pm 0.005$ .

## Procedure

Pipette 2.0 ml of NAD/ATP solution (C), 1.0 ml of Substrate solution (D) and 0.02 ml of G3PDH solution (E) respectively into a quartz cell ( $d = 10\text{mm}$ ) and keep at  $25 \pm 0.5^\circ\text{C}$  for 5 minutes. Then, pipette 0.1 ml of Enzyme solution (G) into the quartz cell and mix well immediately. Keep the reaction mixture at  $25 \pm 0.5^\circ\text{C}$ . Exactly at 2 minutes and 5 minutes after the addition of Enzyme solution (G), measure the absorbances of the reaction mixture at 340 nm ( $A_2$  and  $A_5$ ). As a blank, pipette Diluent (F) into another quartz cell ( $d = 10\text{ mm}$ ) instead of Enzyme solution (G) and take the same procedure described above ( $Ab_2$  and  $Ab_5$ ).

## Calculation

$$\text{Glycerol kinase activity (u/mg)} = \frac{(A_5 - A_2) - (Ab_5 - Ab_2)}{3} \times \frac{1}{6.22} \times 3.12 \times \frac{n}{0.1}$$

- 3: Reaction time
- 6.22: Millimolar absorption coefficient of NADH at 340 nm
- 3.12: Volume of reaction mixture
- 0.1: Volume of Enzyme solution
- n: Dilution factor of Enzyme solution

## Contact

Area	Branch	Location	E-mail
North, Central, South America	Amano Enzyme U.S.A. Co., Ltd.	Illinois, U.S.A.	aeu.sales@amano-enzyme.com
Europe, the Middle East and Africa	Amano Enzyme Europe Limited	Oxfordshire, U.K.	aee.sales@amano-enzyme.com
Asia Pacific	Amano Enzyme Asia Pacific Co., Ltd.	Pathum Thani, THAILAND	aeap.sales@amano-enzyme.com
China	Amano Enzyme Manufacturing (China), Ltd. Shanghai Branch	Shanghai, P.R.CHINA	shanghai@amano-enzyme.com.cn
Japan, Headquarters	Amano Enzyme Inc.	Nagoya, JAPAN	