

Glutaminase "Amano"1

(Glutamine aminohydrolase, EC 3.5.1.2)

Glutaminase "Amano"1 is a glutaminase preparation manufactured by a unique fermentation process using a selected strain of bacteria.

Catalysis

L-Glutamine + H₂O --- L-Glutamic acid + NH₃

Specification and Preparation

Activity: Glutaminase activity Appearance: Light brown powder Additive: Not added

 \geq 800 GTU/g (In-house method)

Characteristics

- 1. Molecular weight: 60,500 (Amino acid sequence) 6.2
- 2. Optimum pH:

3. pH stability: 5.0-9.0 (37°C, 16hrs)

- 4. Optimum temperature: 60-65°C
- 5. Thermal stability: up to 50°C (pH 6.0, 1hr)
- 6. Application: Used for the enzymatic determination of glutamine in serum by coupling with related enzymes in clinical diagnosis.

Expiration (Storage)

18 months from the date of analysis when stored at 20°C or below in a dry place under sealed conditions.

> The information and recommendations contained herein are to the best of our knowledge reliable according to the current scientific and technical level. However, depending upon use method and/or condition, nothing herein is to be construed as a warranty or representation in respect otherwise, including freedom from patent infringement. Users shall make their own test and investigation for their particular purpose. We do not accept any liability for any loss, damage or infringement arising from the use of information and recommendations contained herein.

Safe Handling

- 1. Do not inhale.
- 2. In case of direct contact with skin or eyes, immediately wash or rinse with plenty of water.
- 3. Please refer to SDS for more details.

General properties

The following results demonstrate the activities of enzyme solution prepared in various buffers. Enzyme activity may vary under different experimental conditions.

pH stability pH and Activity Residual activity (%) Relative activity (%) pН pН **Thermostability** Temperature and Activity Residual activity (%) Relative activity (%) With GIn Without GIn Temperature Temperature Stability (powder form) 100.0 Residual activity(%) 80.0 60.0 5°C 40.0 15°C 25°C 20.0



0.0

Months

18 21

Assay method of Glutaminase activity

Principle

L-Glutamine + H₂O $\xrightarrow{\text{Glutaminase}}$ L-Glutamic acid + NH₃ L-Glutamic acid + H₂O + O₂ $\xrightarrow{\text{L-Glutamate oxidase}}$ α -Ketoglutaric acid + NH₃ + H₂O₂ H₂O₂ + TOOS + 4-Aminoantipyrine $\xrightarrow{\text{Peroxidase}}$ Purple pigment

The appearance of purple pigment is measured at 555 nm by spectrophotometry.

Unit Definition

One unit is defined as the enzyme quantity which generates one μ mole of L-glutamic acid per minute under the conditions described below.

Reagents

- A. 1 M Acetic acid buffer (acetic acid sodium acetate, pH6.0)
- B. Triton X-100 solution (50mg/ml)
- C. Sodium hydroxide solution

Weigh 3.0 g of sodium hydroxide in 80 ml of deionized water and cool. Fill up to 100 ml with deionized water.

- D. Perchloric acid solution
 Dilute 8.3 ml of perchloric acid (60%) to 100 ml with deionized water.
- E. Diluent

Mix approx. 600 ml of deionized water, 1 ml of Triton X-100 solution (B), and 10 ml of 1M Acetic acid buffer (A), then fill up to 1000 ml with deionized water.

- F. L-Glutamic acid Assay Kit YAMASA NEO (Yamasa Corporation)
 - F-1: R1 enzyme reagent
 - F-2: R2 enzyme reagent
 - F-3: Glutamic acid standard solution (250 mg/L)
- G. Substrate solution

Weigh 1.00 g of L-glutamine and dissolved in approx. 35ml of deionized water, then add 5 ml of 1 M Acetic acid buffer (A). Fill up to 50 ml with deionized water.

H. Enzyme solution

Weigh out Glutaminase "Amano"1 and dissolve in Diluent (E). Enzyme solution should be prepared so that the value of (As-AsB) becomes in the range of 0.2-0.5 (approx. 0.45-1.10 GTU/ml).



www.amano-enzyme.com

Procedure

<Glutaminase reaction>

Pipette 1.0 ml of Enzyme solution (H) into a test tube and keep at $37 \pm 0.5^{\circ}$ C for 5 minutes. Then, add 1 ml of preheated ($37 \pm 0.5^{\circ}$ C) Substrate solution (G) into the test tube and mix, keep the mixture at $37 \pm 0.5^{\circ}$ C. Exactly after 10 minutes, add Perchloric acid solution (D) and mix, then place in iced water. After 1 minute or more, add 1 ml of Sodium hydroxide solution (C) and mix (Reaction solution).

<L-Glutamic acid Assay Kit>

Pipette 1800 μ I of R1 enzyme reagent (F-1) and 40 μ I of Reaction solution and mix in another test tube, then pipette 1800 μ I of R2 enzyme reagent (F-2) and mix. Exactly after 20 minutes at room temperature, measure the absorbances of the mixture at 555 nm (A_s, Reaction mixture).

<Blank measurement>

As the enzyme blank, pipette 1 ml of Enzyme solution (H) in an empty test tube, add 1 ml of Perchloric acid solution (D) and mix, and keep $37 \pm 0.5^{\circ}$ for 5min. Add 1ml of Substrate solution (G) and immediately place in iced water. After 1 minute or more, add 1ml of Sodium hydroxide solution (C) and mix. After then, follow the procedure of L-Glutamic acid Assay Kit and measure the absorbance at 555nm (A_{SB}).

Also using Glutamic acid standard solution (F-3) or water, measure the absorbance following the same procedure as for the enzyme blank (A_R, A_{RB}).

Calculation

Glutaminase activity (GTU/g) = $\frac{(A_{\rm S} - A_{\rm SB})}{(A_{\rm R} - A_{\rm RB})} \times 250 \times 4 \times \frac{1}{147} \times \frac{1}{10} \times n$

- $A_{S}\!\!:\quad Absorbance \ of \ reaction \ mixture$
- A_{SB}: Absorbance of enzyme blank reaction mixture
- A_R: Absorbance of glutamic acid standard reaction mixture
- A_{RB}: Absorbance of water blank reaction mixture
- 250: Concentration of L-Glutamic acid standard solution
- 4: Volume of reaction mixture
- 147: Molecular weight of L-Glutamic acid
- 10: Reaction time
- N: Dilution factor of Enzyme solution

Contact

Area	Branch	Location	E-mail
North, Central, South America	Amano Enzyme U.S.A. Co., Ltd.	Illinois,U.S.A.	aeu.sales@amano-enzyme.com
Europe, the Middle East and Africa	Amano Enzyme Europe Limited	Oxfordshire, U.K.	aee.sales@amano-enzyme.com
Asia Pacific	Amano Enzyme Asia Pacific Co., Ltd.	Pathum Thani, THAILAND	aeap.sales@amano-enzyme.com
China	Amano Enzyme Manufacturing (China), Ltd. Shanghai Branch	Shanghai, P.R.CHINA	shanghai@amano-enzyme.com.cn
Japan, Headquarters	Amano Enzyme Inc.	Nagoya, JAPAN	

